
Instruction Manual

LIROMA UNIVERSAL ELECTRIC STAND

Max | Multi | Elite | Advanced

Language Index

English: Page 01-22

Nederlands: Pagina 23-42

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Français: Page 63-82

Table of Contents

What's included?	06
Install instructions	08
Panels mounting instructions	10
Panel adjustment instructions	15
Dimensions	18
Control panel instructions	19
Remote control	21
Fault codes and solutions	22

Welcome to Liroma Red Light Therapy

Thank you for choosing Liroma Red Light Therapy.

At Liroma, we understand that your health, recovery, and performance are important. That's why we have developed a range of high-quality red and near-infrared light therapy products designed to support your wellness routine, recovery, and everyday well-being.

Liroma Red Light Therapy products are designed to help support healthy circulation, muscle recovery, joint comfort, relaxation, and overall wellness. Whether you are an athlete, an active individual, or simply looking to support your daily recovery and well-being, red light therapy can be a valuable addition to your self-care routine.

By incorporating red light therapy into your wellness program, you can help support recovery after physical activity, promote relaxation, and contribute to your overall comfort and performance.

Important Notes:

- Liroma Red Light Therapy products do not become excessively hot during normal operation. However, you may experience a gentle warming sensation during treatment. This is completely normal.
- Near-infrared light is invisible to the human eye. Some LEDs may appear inactive during operation, but this is normal and does not affect the performance of the product.

Thank you for trusting Liroma as part of your wellness journey.

With appreciation,

The Liroma Team

WARNINGS:

Please read these instructions carefully before assembling and using the stand.

- Do not assemble or use the stand if any parts or hardware are damaged or missing. If replacement parts are required, please get in touch with Liroma Customer Support for assistance.
- Maximum load capacity: 30 kg (66.1 lbs.).
- Use the stand only as intended by the manufacturer and only with compatible Liroma Red Light Therapy devices.
- This stand contains moving and adjustable components. Use caution during assembly and operation to avoid injury.
- Do not modify the stand or use it in any manner not specified by the manufacturer.
- Liroma is not responsible for damage resulting from unauthorized modifications, improper assembly, misuse, overloading, or damage caused by natural disasters.
- The stand is designed for use on a flat, stable floor surface only. Do not use the stand in any other orientation, such as upside down, on its side, or suspended from the ceiling, as this may cause damage to the stand and the mounted device.
- Keep hands, fingers, hair, and loose clothing away from moving parts during height and angle adjustments.
- Ensure all screws, knobs, and mounting components are securely tightened before installing or operating the light therapy device.


What's Included?

Please check that all items listed below are included before using your Liroma product. If any items are missing or damaged, please get in touch with Liroma customer support before using the product.


Mounting Stand Components & Hardware


A Vertical Support Column x1


B Adjustable Mounting Arm x1


C Base Frame Assembly x1


D Base Cover x1

E Mounting Bracket x2


F Mounting Rail x2

G Locking Castor Wheels x4


H Remote Control x1

I Power Adapter Cable x1

J Cable Clips x4

K Power Extension Cable x1

L Star Knob (Type 1) x2


M Star Knob
(Type 2) x2

N Short Bolts
x8

O Hex Bolts
x4

P Socket Head
Screws x18

Q Flat Head
Screws x10

R Washers
x18


S Allen Key
x1


Instruction
Manual x1


**PrimeFORCE
Max**


**PrimeFORCE
Multi**

**PrimeFORCE
Elite**


**PrimeFORCE
Advanced**


Install Instructions


Step 1 – Install the Wheels.


Step 2 – Install the Support Column


Step 3 – Secure the Support Column


Step 4 – Connect the Cable

Step 5 – Install the Base Cover

Step 6 – Install the Mounting Arm


Panels Mounting Instructions

1. PrimeFORCE Max Mounting

Step 1:

Attach the two mounting rails (F x2). Secure with screws.


Step 2:

Mount onto the mounting arm (B).

Step 3:


Secure with the mounting screws.

2. PrimeFORCE Multi Mounting


Step 1:

Remove the feet from one panel and place both panels face down on a flat surface.


Step 2:

Attach the mounting brackets (F x2) with supplied screws.

Step 3:


Install on the mounting arm (B). Tighten the screws.

3. PrimeFORCE Elite Mounting


Step 1:

Attach the mounting rails (F x2). Secure with screws.

Step 2:


Secure with the mounting screws.

4. PrimeFORCE Advanced Mounting


Step 1:

Attach the mounting bracket (E) using the supplied screws.

Step 2:


Secure with the mounting screws.


Panel Adjustment Instructions

Angular Adjustment


Note: Only when the bulb is in a downward-tilted position, follow the instructions on the next page to adjust it to a horizontal position. If the bulb is already in a horizontal state, no adjustment is necessary.


1. Loosen 4 screws by 2–3 millimeters.

2. Raise the lamp up, insert the gasket between the two screws near the lamp.


3. Tighten all the screws.


Panel Rotation Adjustment


1. Raise the stand to its highest position.
2. Hold the panel with one hand, raise the handle upwards.
3. Rotate the panel to your desired position.

4. Press the handle down to lock it.


Dimensions

Control Panel Instructions

Control Panel

Initial Power-On

1. Connect the stand to a 100–240 V AC power supply. The buzzer will emit a single beep, and the controller will display the current height.
2. If the height value shown on the controller does not match the actual height, or if the height has changed after reassembly, a reset operation will be necessary. Otherwise, the system can be used directly with the factory settings.

Reset Operation

A reset operation is required after reassembly or whenever the displayed height is incorrect.

1. Press and hold the ▼ button. After approximately 5 seconds, the display will show "rst".
2. Continue holding the ▼ button while the stand slowly descends to its lowest point.
3. Once the stand reaches the lowest position and rebounds approximately 1 cm (0.4 in), release the ▼ button to conclude the reset operation.
4. If an emergency arises during the reset, you can release the ▼ button at any.
5. time to stop the stand and end the reset.
6. After an interrupted reset, it is necessary to reinitiate the reset operation until it is successfully completed; no other operations can be performed until then.

Important: If the reset process is interrupted, it must be restarted and completed successfully before any other functions can be used.

Height Adjustment

- **Raising the Stand:** Press and hold the ▲ button to raise the height. Release the button to stop.
- **Lowering the Stand:** Press and hold the ▼ button to lower the height. Release the button to stop.

Preset Height

1. Press the [S] button, [S] appears on the control display.
2. Press [1] or [2] or [3] button to save the height to the selected button.
3. The control can store up to three heights, preserved even after power outages.
4. To recall a saved height, press the corresponding button. The stand will move automatically to that height.
5. Press any button during movement to stop immediately.

Lock / Unlock


1. Press and hold the ▲ (Up) and ▼ (Down) buttons simultaneously for 5 seconds.
2. The display will show "LOC", the stand is locked and cannot be adjusted.
3. Press and hold the same buttons until the display returns to the normal digital readout, the stand is unlocked.

Energy Saving Mode

If no buttons are pressed for 1 minute, the stand enters Energy Saving Mode. The display will turn off to reduce power consumption.
Press any button to exit this mode.

Remote Control

The remote control is pre-programmed and ready for use. Pairing is NOT NEEDED unless you have to replace the remote.

Reset

1. Press and hold the **RST** button for 3 seconds a beep sound indicates that the reset process has started.
2. The stand lowers to the minimum height and rebounds 1cm.
3. Two beeps sound to signal reset completion.

Pairing

Pairing is only required when replacing the remote control.

1. Unplug the power cable of the stand.
2. Press and hold the **SET** button on the Remote Controller without releasing it.
3. Power on the stand, and upon hearing a "beep" sound indicating that pairing is in progress, release the **SET** button.

▲ Press and hold to raise the height, release to stop.

▼ Press and hold to lower the height, release to stop.

M1 Memorized height 1

M2 Memorized height 2

Preset Height

1. The stand stops at an ideal height.
2. Press and hold **M1** or **M2** button for 3 seconds until a "beep" sound. To save the current height.
3. Press **M1** or **M2** button to the memorized height.

Fault Codes and Solutions

Fault	Codes	Solutions
Overheat Protection	HOT	<p>Solution 1: No action needed. Automatically recover in 18 minutes.</p> <p>Solution 2: Disconnect the power, wait for the screen to turn off, then reconnect the power and resume use.</p>
Reset Operation Not Completed	rst	Press and hold the ▲ button on the control panel until the reset is fully completed.
Overload Protection	E11	Remove any heavy objects from the stand. Wait for the alarm to clear before using it again.
Motor Operation Abnormality	E12	Turn off the power, check the connection between the motor and the controller. Turn it on again.
Anti-Collision Protection	E14	It will automatically reset by running in reverse for a specific distance. It should operate normally after removing any obstacles.
Motor Stall Alarm	E15	If there is a severe overload or the transmission mechanism is stuck, remove any overloaded objects and check for foreign objects within the transmission mechanism.
Motor Wire Not Connected	E16	Check and reconnect the connection line between the motor and the controller.

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